

# Which FOTODYNE Imaging System is Right for Your Lab?

## FOTO/Analyst®

### Apprentice

These entry-level systems combine high resolution color digital cameras and inexpensive photographic hoods. Auto Connect USB hook up makes them compatible with PCs and Macs.

\$1,695 - \$3,490



## FOTO/Analyst® Express

Live images and greater sensitivity provided by CCD-based cameras make these hood-mount systems a step up from their Apprentice counterparts. Combine with your computer to acquire, print and store digital images. Available in PC and Mac formats. Great for labs that want to break away from a shared system. Network compatible. Several camera options, including color.

\$3,400 - \$6,950



## FOTO/Analyst®

### MiniVisionary

These stand-alone introductory-level systems are designed to simply create quick, inexpensive thermal prints. No computer required. Ease-of-use and rapid throughput make these systems popular in large teaching labs and QC labs.

\$4,200 - \$5,200



## FOTO/Analyst® Investigator

Our most popular system. PC or Mac compatible workstations with benchtop darkroom and transilluminator. Acquire, print and store images of a variety of sample types and sizes. Ideal for small research labs with dedicated users and for labs that wish to use their imaging system on a network. Several camera options, including color.

\$4,800 - \$10,650



## FOTO/Analyst® Investigator/FX

The flexible FOTO/FX benchtop darkroom makes this our most customizable product line. Choose system features based on your needs - camera specifications, manual or motorized zoom lens, transilluminator size, illumination wavelengths, and more. The Investigator/FX can easily be upgraded with a cooled camera for chemiluminescent samples. PC, I-Mac, and network compatible.

\$8,500 - \$14,250



## FOTO/Analyst® Luminary FX

Top-of-the-line cooled CCD camera workstations for obtaining digital images of chemiluminescent samples, gels, petri plates, microassay plates and more. Each workstation features the flexible FOTO/FX benchtop darkroom, transilluminator and filters. PC, I-Mac, and network compatible. Three camera choices.

\$15,100 - \$32,500

Call or e-mail today to receive more information, request a price quotation, or schedule a demo.  
customer.service@fotodyne.com • 1-800-DNA-FOTO (362-3686) • www.fotodyne.com  
FOTODYNE Incorporated • 950 Walnut Ridge Drive • Hartland, WI 53029-9388 • USA

**FOTODYNE**  
INCORPORATED